

## Protocol for Sequencing BAC subclones

### 1. Day 1

- a. Prepare Media and agar plates needed for the next week of steps.
  - i. 1 sleeve of LB agar plates with Chloramphenicol (12.5mg/ml), X-Gal (20mg/ml) and IPTG (200mg/ml). May only need 1 plate depending on number of BACs being cloned.
  - ii. 2 sleeves of LB agar plates with Kanamycin (10mg/ml), X-Gal (20mg/ml) and IPTG (200mg/ml).
  - iii. 5ml LB media tubes (1 per BAC)
  - iv. 200ml LB media flasks (1 per BAC)
  - v. Solutions not included in Qiagen Large Construct Kit
  - vi. Solutions not included in Invitrogen Topo Shotgun Subcloning Kit
  - vii. Woo's Freezing Broth
  - viii. Terrific Broth
  - ix. Solution I, Solution II and Solution III (prepare morning of use)

### 2. Day 2

- a. Pull 1ul of needed BAC from library plate and put into 50-100ul ddH<sub>2</sub>O.
- b. Streak 10-15ul of BAC onto LB agar plates with Chloramphenicol (12.5mg/ml), X-Gal (20mg/ml) and IPTG (200mg/ml).
- c. Let plate grow at 37°C overnight (approximately 18 hours) and then place in 4°C for a couple hours to allow for blue/white selection.

### 3. Day 3

- a. In the morning, start a 5ml culture using 5ml LB media with 5ul Chloramphenicol. Pick a single white colony off the plate, using a pipette tip, and eject the tip into the LB media.
- b. Allow this to grow for 6-8 hours at 37°C 280rpm.
- c. Transfer 500ul of 5ml culture into 100-500ml LB media with 100-500ul Chloramphenicol. Allow this to grow for 18 hours at 37°C 280rpm.

### 4. Day 4

- a. Perform DNA extraction using the Qiagen Large Construct Kit, following the kit protocol.
- b. Determine DNA concentration via the nanodrop.
- c. May also want to determine DNA concentration by digesting the BAC with NotI and running the digest out on a CHEF gel.

### 5. Day 5

- a. Prepare the DNA for subcloning following the Invitrogen TOPO Shotgun Subcloning Kit protocol. \*\*
  - i. Shear DNA to approximately 3kb fragments.
  - ii. Purify Sheared DNA (not needed but may be helpful)
  - iii. Blunt-End repair DNA
  - iv. Dephosphorylate DNA
  - v. TOPO cloning of DNA
  - vi. Spread cloned cells onto LB agar plates with Kanamycin (10mg/ml), X-Gal (20mg/ml) and IPTG (200mg/ml). I spread 75ul of cells onto 1 plate and let grow overnight to make sure transformation worked, then the next day spread 100ul onto 4 plates.
  - vii. Let plates dry in hood and then place in 37°C incubator for 18 hours.

\*\*These steps may take more than 1 day, I usually Shear, Purify, Blunt-End Repair and start Dephosphorylation on the first day (to step 4) then finish Dephosphorylation and clone on the second day.

6. Day 6
  - a. Prepare 384-well plates for colonies by adding 70ul Woo's Freezing Broth with Kanamycin (10mg/ml) to each well. If using a lot of plates, use Q-Fill to load.
  - b. Pick white colonies from the LB agar plates containing subcloned DNA into the 384-well plate, using sterile toothpicks. Put one colony per well.
  - c. Allow plate to grow for 18-22 hours at 37°C.
7. Day 7
  - a. Fill 4-96 deep well blocks with 300ul Terrific Broth and Kanamycin (10mg/ml) into each well.
  - b. Replicate the 384-well plate into the blocks using the 96 long-pin replicator.
    - i. Fill the sand with 95% Ethanol.
    - ii. Push replicator into sand a couple times to help clean the pins.
    - iii. Flame replicator with ethanol burner and allow to cool on bench.
    - iv. Place replicator into block to further cool then transfer to 384-well plate placing the upper left pin into well A1.
    - v. Swish replicator in plate a few times to collect cells and return to 96-well block and swish to release cells.
    - vi. Push replicator into the sand and repeat steps ii-v for all 4 blocks replicating the 384-well plate into the 4-96 well blocks (A1, A2, B1, B2).
  - c. Cover the blocks with 2 layers of AeraSeal and place in the 37°C shaking incubator for 18 hours at 280rpm.
  - d. Cover the 384-well plate with Foil Seal and place into -80°C for future use.
8. Day 8
  - a. Perform 96 well High Copy Plasmid Preps
    - i. Pellet the cells by centrifuging in refrigerated bench-top centrifuge at 3000rpm for 15 min, 4°C.
    - ii. While centrifuging, label 1 Uniplate and 1 Unifilter for each block.
    - iii. Pour off supernatant and place blocks upside to dry on paper towels.
    - iv. Add 150ul of cold Solution I to each well, Cover with tape and vortex to re-suspend pellets then remove tape. The culture should only be in Solution I for 10 minutes.
    - v. Add 150ul of Solution II to each well. Cover with tape and mix by inverting 3-4 times then remove tape. The culture should only be in Solution II for 5 minutes.
    - vi. Add 150ul of cold Solution III to each well. Cover with tape and mix by inverting 3-4 times.
    - vii. Place block on ice for at least 10 minutes.
    - viii. Add 270ul of isopropanol (0.6 volumes) to each well of the Uniplate.
    - ix. Place Unifilter on Uniplate and add 450ul of each sample into the Unifilter.  
\*\*\*Make sure blocks are all in the same direction\*\*\*
    - x. Centrifuge at 3200 rpm for 30 minutes 4°C.
    - xi. Throw out Unifilters, pour off supernatant and let drain upside down.
    - xii. Add 270ul of 70% Ethanol to each well.
    - xiii. Centrifuge at 3200 rpm for 10 minutes 4°C.
    - xiv. Pour off supernatant and allow to drain upside down. Place Uniplates in hood to dry completely for approximately 10 minutes.
    - xv. Add 100ul ddH<sub>2</sub>O to each well and shake Uniplate slightly to resuspend pellet.
    - xvi. Cover with Thermoseal and place in -20C until further use.
9. Day 9
  - a. Set-up cycle sequencing reactions using ABI Big Dye v3.1 protocol.
    - i. Allow Subclone samples to thaw in Uniplate. Label half-skirted PCR plates with BAC name, Uniplate name and primer for sequencing.

- ii. Add 1ul sample from Uniplate into half-skirted PCR plate.
- iii. Mix cocktail of sequencing reagents (10ul total volume per well). Make 100 reactions per plate.
  - 1. 3ul ddH<sub>2</sub>O 300ul in cocktail
  - 2. 3ul 5x Sequencing Buffer 300ul
  - 3. 2ul Forward or Reverse M13 primer (1.6pmol/ul) 200ul
  - 4. 1ul ABI Big Dye Ready Reaction Mix 100ul
- iv. Aliquot 9ul per well, tap plate to mix then quick centrifuge to make sure all contents are at the bottom of plate.
- v. Place plate into Thermal Cycler and run on subclone sequence program.
  - 1. Denature at 96°C for 1 minute
  - 2. Denature at 96°C for 10 seconds
  - 3. Anneal at 50°C for 5 seconds
  - 4. Extend at 60°C for 4 minutes
  - 5. Repeat steps 2-4 24 times
  - 6. Place at 4°C forever, until plates are removed.

#### 10. Day 10

- a. Purify cycle sequencing using a EDTA/Ethanol precipitation protocol.
  - i. Add 5ul of 125mM EDTA to each well. Tap plate to mix.
  - ii. Add 60ul of 100% Ethanol to each well. Cover with seal and invert plate to mix.
  - iii. Let plate sit at Room Temperature for 15 minutes to precipitate DNA.
  - iv. Centrifuge plate at 3200xg for 30 minutes, 4°C.
  - v. Remove plates and discard supernatant.
  - vi. Add 60ul of 70% Ethanol to each well.
  - vii. Centrifuge at 1650xg for 15 minutes, 4°C.
  - viii. Remove plates, discard supernatant and invert plates on paper towels.
  - ix. Centrifuge inverted plates at 185xg for 1 min to remove ethanol.
  - x. Dry samples in hood for 10 minutes then resuspend pellets with 10ul Formamide.
  - xi. Cover with Septa Seal and place into Sequencer holders. If these are not available, wrap plates in Aluminum foil and place in 4°C.
- b. Enter sequence input into Excel Spreadsheet
- c. Put the Sequence input files into Greg's personal file folder→Sequencer3730→Sequencer input→Lori→create folder with date created.
- d. Give plates to Greg to sequence.