

DNA EXTRACTIONS FOR PCR

Extraction buffer-50 mls

| | | |
|---------------------|-------------|--------|
| 2M Sorbitol | 3.5 mls | 0.14M |
| 1M Tris pH 7.5 | 11.0 mls | 0.22M |
| 0.5M EDTA pH 8.0 | 2.2 mls | 0.22M |
| 5M NaCl | 8.0 mls | 0.8M |
| CTAB | 0.4 gms | 0.022M |
| N-lauroyl. sarcosyl | 0.5 gms | 0.034M |
| dd H ₂ O | 25 mls | |
| Na bisulfate | 0.1g/30 mls | |

Where volumes are specified assume that 1 ml is equal to 1 gram. Weigh solutions into a new sterile tissue culture bottle. It is best done the night before so that the CTAB and sarcosyl have plenty of time to dissolve. Usually I leave it in the tissue culture growth rooms. Approximately 2 mls of buffer are required for each sample.

1. Collect sample into plastic bag (1 leaf tip for onion).
2. Roll sample with a wall paper roller or lead pipe and then add 1.5 mls of extraction buffer (2X 750 μ ls). Roll again with the roller to ensure the buffer and plant material are well mixed.
3. Pipette 300 μ l of chloroform into tube. Transfer sample from the bag into the tube by pouring from one corner.
4. Invert to mix and incubate at 65°C for 30 minutes.
5. Invert tube to ensure chloroform and sample are well mixed. Centrifuge at 13,000 rpm for 5 minutes.
6. Pipette 600 μ l of isopropanol into a new tube. Pipette up to 1 ml of supernatant into this tube.
7. Invert to mix and centrifuge at 13,000 rpm for 5 minutes to pellet DNA.
8. Pour off supernatant and wash pellet with 1 ml of 70% ethanol (room temperature).
9. Optional centrifugation after the 70% wash if the pellet has become dislodged; 13,000 rpm for 2-3 minutes.
10. Dry pellet in laminar flow unit. Tubes dry faster if laid on their side.
11. When pellets are nice and clear, resuspend in 200-400 μ l PCR ddH₂O. (50 -100 μ TE for onion?)
12. Quantitate DNA on a 1% agarose gel. Dilute 1 μ l of sample into a mixture of 15 μ l dd H₂O + 4 μ l of sample buffer. Mix well and load 5 μ l on the gel (= 1/4 of a μ l). NB: Samples must be completely dissolved before quantitation.